

Aditamentos

DIO IMPLANT


Narrow

Diámetro: Ø3.0/Ø3.3

Longitud: 8.5/10/11.5/13/15

Ø3.0 Accesorio: Solo para anterior

Superficie del Implante: HSA/UV Active



Cement-Retained Restorations
Temporary / Cemented / Angled Abutment

Screw-Retained Restorations
UCLA CCM Abutment

Cement-Retained Restorations
Solid Abutment




Stainless Steel Housing




OT Cap




Ball Abutment Analog




Ball Abutment

Overdenture-Retained Restorations
Ball Abutment


Cylinder Screw




Gold Cylinder


Plastic Cylinder


Temporary Cylinder




Abutment Analog


**Pick-up
Impression Coping**


**Transfer
Impression Coping**


**Healing
Cap**


**Conical
Abutment**


**Multiunit
Abutment**

Screw-Retained Restorations
Multiunit / Conical Abutment


UF II Narrow Fixture


Cover Screw


Healing Abutment

Narrow Fixture

Hex	Length	Ø3.0	Ø3.3
1.7	8.5	UF(II)N 3008S	UF(II)N 3308S
	10	UF(II)N 3010S	UF(II)N 3310S
	11.5	UF(II)N 3011S	UF(II)N 3311S
	13	UF(II)N 3013S	UF(II)N 3313S
	15	UF(II)N 3015S	UF(II)N 3315S
Apex		Ø1.6	Ø1.9

Healing Abutment

Height	Cuff	Ø4.1	Ø4.2	Ø4.6	Ø4.7
2	0		UNSHA 4020		UNSHA 4520
4	2		UNSHA 4024		UNSHA 4524
5.5	3		UNSHA 4035		UNSHA 4535
7	4		UNSHA 4047		UNSHA 4547
8.5	5		UNSHA 4058		UNSHA 4558
10	6		UNSHA 4060		UNSHA 4560

Temporary Abutment

Cuff	Ø4.0	
	Hex (Guide Pin)	Non-Hex
1	UNSTA 4010H	UNSTA 4010N
3	UNSTA 4030H	UNSTA 4030N

Cemented Abutment

Cuff	Ø4.0			
	Hex		Non-Hex	
	L 5.5	L 7	L 5.5	L 7
1	UNSCA 40105H	UNSCA 40107H	UNSCA 40105N	UNSCA 40107N
2	UNSCA 40205H	UNSCA 40207H	UNSCA 40205N	UNSCA 40207N
3	UNSCA 40305H	UNSCA 40307H	UNSCA 40305N	UNSCA 40307N
4	UNSCA 40405H	UNSCA 40407H	UNSCA 40405N	UNSCA 40407N
5	UNSCA 40505H	UNSCA 40507H	UNSCA 40505N	UNSCA 40507N
6	UNSCA 40605H	UNSCA 40607H	UNSCA 40605N	UNSCA 40607N
7	UNSCA 40705H	UNSCA 40707H	UNSCA 40705N	UNSCA 40707N

Angled Abutment

Cuff	Ø4.0	
	Hex B	Non-Hex
2	UNSA A 402015BH	UNSA A 402015N
4	UNSA A 404015BH	UNSA A 404015N

UCLA CCM Abutment

Cuff	Ø4.0	
	Hex	Non Hex
1	UNSCCA 4010H[H]	UNSCCA 4010N[H]
3	UNSCCA 4030H[H]	UNSCCA 4030N[H]

Solid Abutment

Cuff	Ø4.0	
	L 5.5	L 7
1	UNSSA 40105	UNSSA 40107
2	UNSSA 40205	UNSSA 40207
3	UNSSA 40305	UNSSA 40307
4	UNSSA 40405	UNSSA 40407
5	UNSSA 40505	UNSSA 40507
6	UNSSA 40605	UNSSA 40607
7	UNSSA 40705	UNSSA 40707
Protect Cap	UNSAAC 4005	UNSAAC 4007

Ball Abutment

Length	Ø4.0
1	UNBSA 3510
2	UNBSA 3520
3	UNBSA 3530
4	UNBSA 3540
5	UNBSA 3550
6	UNBSA 3560
7	UNBSA 3570

Multiunit Abutment

Hex

Cuff	Angle 20°
2.5	UNSMUA 482520H
3	UNSMUA 483020H
4	UNSMUA 484020H
5	UNSMUA 485020H

Non-Hex

Cuff	Angle 20°
2.5	-
3	UNSMUA 483020N
4	UNSMUA 484020N
5	UNSMUA 485020N

Regular

Wide

Diámetro: Ø3.8 / Ø4.0 / Ø4.5 / Ø5.0 / Ø5.5
Longitud: 7 / 8.5 / 10 / 11.5 / 13
Superficie del Implante: HSA / UV Active / MG

Diámetro: Ø6.0 / Ø6.5 / Ø7.0
Longitud: 7 / 8.5 / 10 / 11.5
Superficie del Implante: HSA / UV Active / MG+



Abutment Screw



Temporary Abutment



Cemented Abutment



Mill Abutment



Angled Abutment



UCLA Gold/CCM Abutment



UCLA Plastic Abutment



Fixture Analog



Pick-up Impression Coping



Transfer Impression Coping

Cement-Retained Restorations
Temporary / Cemented / Mill / Angled Abutment

Screw-Retained Restorations
UCLA Gold / UCLA CCM / UCLA Plastic Abutment



UF II
Regular Fixture



UF II
Wide Fixture



Cover Screw



Healing Abutment



O-Ring



Stainless Steel Housing



Denture Cap



Retainer



Ball Cap



OT Cap



Dalbo Plus



Male



Ball Abutment Analog



Ball Abutment



Locator

Overdenture-Retained Restorations
Ball Abutment / Locator



UF II Regular Fixture



UF II Wide Fixture



Cover Screw




Healing Abutment


Cylinder Screw



Octa Abutment Analog




Screw-Retained Restorations
Octa Abutment



Cylinder Screw



Abutment Analog



Screw-Retained Restorations
Multiunit / Conical Abutment



Regular Fixture

HSA		Surface: HSA, UV Active, MG ⁺				
Hex	Length	Ø3.8	Ø4.0	Ø4.5	Ø5.0	Ø5.5
2.5	7	-	-	UF(II) 45075	UF(II) 50075	UF(II) 55075
	8.5	UF(II) 38085	UF(II) 40085	UF(II) 45085	UF(II) 50085	UF(II) 55085
	10	UF(II) 38105	UF(II) 40105	UF(II) 45105	UF(II) 50105	UF(II) 55105
	11.5	UF(II) 38115	UF(II) 40115	UF(II) 45115	UF(II) 50115	UF(II) 55115
	13	UF(II) 38135	UF(II) 40135	UF(II) 45135	UF(II) 50135	UF(II) 55135
	15	UF(II) 38155	UF(II) 40155	UF(II) 45155	UF(II) 50155	UF(II) 55155
	16	UF(II) 38165	UF(II) 40165	UF(II) 45165	UF(II) 50165	UF(II) 55165
18	UF(II) 38185	UF(II) 40185	UF(II) 45185	UF(II) 50185	UF(II) 55185	
Apex		Ø2.1	Ø2.2	Ø2.6	Ø3.1	Ø3.7

Wide Fixture

HSA		Surface: HSA, UV Active		
Hex	Length	Ø6.0	Ø6.5	Ø7.0
2.5	7	UF(II) 60075	UF(II) 65075	UF(II) 70075
	8.5	UF(II) 60085	UF(II) 65085	UF(II) 70085
	10	UF(II) 60105	UF(II) 65105	UF(II) 70105
	11.5	UF(II) 60115	UF(II) 65115	UF(II) 70115
	13	UF(II) 60135	UF(II) 65135	UF(II) 70135
15	UF(II) 60155	UF(II) 65155	UF(II) 70155	
Apex		Ø5.4	Ø5.9	Ø6.4

Healing Abutment

Cuff		Ø4.0			
2		SSHA 4020			
4		SSHA 4040			
6		SSHA 4060			
8		SSHA 4080			
Height	Cuff	D1 - Ø4.6	D1 - Ø5.6	D1 - Ø6.6	D1 - Ø7.6
		D2 - Ø4.5	D2 - Ø5.5	D2 - Ø6.5	D2 - Ø7.5
2	2	SSHA 4520	SSHA 5520	SSHA 6520	SSHA 7520
4	2	SSHA 4524	SSHA 5524	SSHA 6524	SSHA 7524
5.5	3	SSHA 4535	SSHA 5535	SSHA 6535	SSHA 7535
7	4	SSHA 4547	SSHA 5547	SSHA 6547	SSHA 7547
Height	Cuff	D1 - Ø8.2		D1 - Ø9.2	
		D2 - Ø8.0		D2 - Ø9.0	
2	1	SSHA 8212		SSHA 9212	
3	2	SSHA 8223		SSHA 9223	
4	3	SSHA 8234		SSHA 9234	
5	4	SSHA 8245		SSHA 9245	

Temporary Abutment

Cuff	Ø4.5	
	Hex	Non-Hex
1	SSTA 4510H[H]	SSTA 4510N[H]
3	SSTA 4530H[H]	SSTA 4530N[H]

Cemented Abutment

Hex Type

Length	Cuff	Ø4.5	Ø5.5	Ø6.5	Ø7.5
4	1	SSCA 45104H(II)[H]	SSCA 55104H(II)[H]	SSCA 65104H(II)[H]	SSCA 75104H(II)[H]
	2	SSCA 45204H(II)[H]	SSCA 55204H(II)[H]	SSCA 65204H(II)[H]	SSCA 75204H(II)[H]
	3	SSCA 45304H(II)[H]	SSCA 55304H(II)[H]	SSCA 65304H(II)[H]	SSCA 75304H(II)[H]
	4	SSCA 45404H(II)[H]	SSCA 55404H(II)[H]	SSCA 65404H(II)[H]	SSCA 75404H(II)[H]
5.5	5	SSCA 45504H(II)[H]	SSCA 55504H(II)[H]	SSCA 65504H(II)[H]	SSCA 75504H(II)[H]
	1	SSCA 45105H(II)[H]	SSCA 55105H(II)[H]	SSCA 65105H(II)[H]	SSCA 75105H(II)[H]
	2	SSCA 45205H(II)[H]	SSCA 55205H(II)[H]	SSCA 65205H(II)[H]	SSCA 75205H(II)[H]
	3	SSCA 45305H(II)[H]	SSCA 55305H(II)[H]	SSCA 65305H(II)[H]	SSCA 75305H(II)[H]
	4	SSCA 45405H(II)[H]	SSCA 55405H(II)[H]	SSCA 65405H(II)[H]	SSCA 75405H(II)[H]
7	5	SSCA 45505H(II)[H]	SSCA 55505H(II)[H]	SSCA 65505H(II)[H]	SSCA 75505H(II)[H]
	1	SSCA 45107H(II)[H]	SSCA 55107H(II)[H]	SSCA 65107H(II)[H]	SSCA 75107H(II)[H]
	2	SSCA 45207H(II)[H]	SSCA 55207H(II)[H]	SSCA 65207H(II)[H]	SSCA 75207H(II)[H]
	3	SSCA 45307H(II)[H]	SSCA 55307H(II)[H]	SSCA 65307H(II)[H]	SSCA 75307H(II)[H]
	4	SSCA 45407H(II)[H]	SSCA 55407H(II)[H]	SSCA 65407H(II)[H]	SSCA 75407H(II)[H]
Screw		SSCA 45507H(II)[H]	SSCA 55507H(II)[H]	SSCA 65507H(II)[H]	SSCA 75507H(II)[H]
		SSC 2008SH	SSC 2008H	SSC 2008H	SSC 2008H

Non-Hex Type

Length	Cuff	Ø4.5	Ø5.5	Ø6.5	Ø7.5
4	1	SSCA 45104N(II)[H]	SSCA 55104N(II)[H]	SSCA 65104N(II)[H]	SSCA 75104N(II)[H]
	2	SSCA 45204N(II)[H]	SSCA 55204N(II)[H]	SSCA 65204N(II)[H]	SSCA 75204N(II)[H]
	3	SSCA 45304N(II)[H]	SSCA 55304N(II)[H]	SSCA 65304N(II)[H]	SSCA 75304N(II)[H]
	4	SSCA 45404N(II)[H]	SSCA 55404N(II)[H]	SSCA 65404N(II)[H]	SSCA 75404N(II)[H]
	5	SSCA 45504N(II)[H]	SSCA 55504N(II)[H]	SSCA 65504N(II)[H]	SSCA 75504N(II)[H]
5.5	1	SSCA 45105N(II)[H]	SSCA 55105N(II)[H]	SSCA 65105N(II)[H]	SSCA 75105N(II)[H]
	2	SSCA 45205N(II)[H]	SSCA 55205N(II)[H]	SSCA 65205N(II)[H]	SSCA 75205N(II)[H]
	3	SSCA 45305N(II)[H]	SSCA 55305N(II)[H]	SSCA 65305N(II)[H]	SSCA 75305N(II)[H]
	4	SSCA 45405N(II)[H]	SSCA 55405N(II)[H]	SSCA 65405N(II)[H]	SSCA 75405N(II)[H]
	5	SSCA 45505N(II)[H]	SSCA 55505N(II)[H]	SSCA 65505N(II)[H]	SSCA 75505N(II)[H]
7	1	SSCA 45107N(II)[H]	SSCA 55107N(II)[H]	SSCA 65107N(II)[H]	SSCA 75107N(II)[H]
	2	SSCA 45207N(II)[H]	SSCA 55207N(II)[H]	SSCA 65207N(II)[H]	SSCA 75207N(II)[H]
	3	SSCA 45307N(II)[H]	SSCA 55307N(II)[H]	SSCA 65307N(II)[H]	SSCA 75307N(II)[H]
	4	SSCA 45407N(II)[H]	SSCA 55407N(II)[H]	SSCA 65407N(II)[H]	SSCA 75407N(II)[H]
	5	SSCA 45507N(II)[H]	SSCA 55507N(II)[H]	SSCA 65507N(II)[H]	SSCA 75507N(II)[H]
Screw		SSC 2008SH	SSC 2008H	SSC 2008H	SSC 2008H

CCM Type

Cuff	Ø4.0		Ø4.5	
	Hex	Non Hex	Hex	Non Hex
1	CCSS 4010H[H]	CCSS 4010N[H]	CCSS 4510H[H]	CCSS 4510N[H]
3			CCSS 4530H[H]	CCSS 4530N[H]

Solid Abutment

Length	Cuff	Ø4.5	Ø5.5	Ø6.5	Ø7.5
4	1	SSSA 45104	SSSA 55104	SSSA 65104	SSSA 75104
	2	SSSA 45204	SSSA 55204	SSSA 65204	SSSA 75204
	3	SSSA 45304	SSSA 55304	SSSA 65304	SSSA 75304
	4	SSSA 45404	SSSA 55404	SSSA 65404	SSSA 75404
	5	SSSA 45504	SSSA 55504	SSSA 65504	SSSA 75504
5.5	1	SSSA 45105	SSSA 55105	SSSA 65105	SSSA 75105
	2	SSSA 45205	SSSA 55205	SSSA 65205	SSSA 75205
	3	SSSA 45305	SSSA 55305	SSSA 65305	SSSA 75305
	4	SSSA 45405	SSSA 55405	SSSA 65405	SSSA 75405
	5	SSSA 45505	SSSA 55505	SSSA 65505	SSSA 75505
7	1	SSSA 45107	SSSA 55107	SSSA 65107	SSSA 75107
	2	SSSA 45207	SSSA 55207	SSSA 65207	SSSA 75207
	3	SSSA 45307	SSSA 55307	SSSA 65307	SSSA 75307
	4	SSSA 45407	SSSA 55407	SSSA 65407	SSSA 75407
	5	SSSA 45507	SSSA 55507	SSSA 65507	SSSA 75507

Ball Abutment

Length	Ø3.5
1	SSBA 3510
2	SSBA 3520
3	SSBA 3530
4	SSBA 3540
5	SSBA 3550
6	SSBA 3560

Multiunit Abutment

Hex


Cuff	Angle 20°	Angle 30°
2.5	UMA 482520H	-
3	UMA 483020H	UMA 483030H
4	UMA 484020H	UMA 484030H
5	UMA 485020H	UMA 485030H

Non-Hex

Cuff	Angle 20°	Angle 30°
2.5	UMA 482520N	-
3	UMA 483020N	UMA 483030N
4	UMA 484020N	UMA 484030N
5	UMA 485020N	UMA 485030N

KIT

UF II Int. Submerged



Torque Wrench
DTW 0060



Open Wrench
OW 004



Depth Gauge
DPG 2050



UF II Wide Kit



UF II Narrow Kit



UF II Master Kit



Sinus Master Kit



Dr.SOS+ Kit

Aditamentos